CryoSure-DEX40

- (EN)** Cryoprotectant for the cryopreservation of hematopoietic stem cells for transplantation
- (DE)** Kryoprotektivum zur Kryokonservierung hämatopoetischer Stammzellen, die zur Transplantation bestimmt sind
- (F)** Solution cryoprotectrice employée pour la cryoconservation des cellules souches hématopoïétiques à visée thérapeutique

(EN) Sterile acc.to EP/USP	Pyrogen-free acc. to EP	Endotoxin-free acc. to EP/USP	Free of mycoplasma acc. to EP
(DE) Steril gem. EP/USP	Pyrogenfrei gem. EP	Endotoxinfrei gem. EP/USP	Frei von Mykoplasmen gem. EP
(F) Stérile selon Ph. Eur./USP	Exempt de pyrogènes selon Ph. Eur.	Exempt d'endotoxines bactériennes selon Ph. Eur./USP	Exempt de mycoplasmes selon Ph. Eur.

(EN) Ingredients:	55 g/dl DMSO USP Grade, 5 g/dl Dextran 40 USP Grade, water for injection EP
(DE) Inhaltsstoffe:	55 g/dl DMSO USP Grade, 5 g/dl Dextran 40 USP Grade, Wasser für Injektionszwecke EP
(F) Ingrédients:	55 g/dl DMSO USP Grade, 5 g/dl Dextrane 40 USP Grade, eau pour préparations injectables selon Ph. Eur

CE
0482

WARNING: Do not autoclave. **Do not use unless solution is clear.** **Not for injection.**
ACHTUNG: Nicht autoklavieren. **Nur klare Lösung verwenden.** **Nicht injizieren.**
ATTENTION: Ne pas passer en autoclave. **Seulement utiliser la solution claire.** **Ne pas injecter.**

STERILE A

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Introduction:

- (EN) CryoSure-DEX40 is employed in the course of longterm storage and cell counting of hematopoietic stem cell preparations.

CryoSure-DEX40 is a ready-to-use cryoprotective solution for the addition to a buffy coat suspension from volume-reduced cord blood according to the method of Rubinstein et al (1). CryoSure-DEX40 is a solution consisting of 50% v/v DMSO and 50% v/v of a 10% aqueous solution of Dextran 40.

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EN It contains the cryoprotectants DMSO and Dextran 40, which have their effects inside and outside the cell, respectively.

DMSO (Dimethyl Sulfoxide) is a cryoprotectant which penetrates the cell wall and takes its cryoprotectant effect within the cell. It reduces the osmotic stress on the cells during freezing and thawing (2, 3, 4, 5, 6) and antagonizes the osmotic shock (7). Also DMSO protects the cells by reducing dehydration and shrinkage of the cells during the freezing process (5, 8). After thawing DMSO has to be removed from the stem cell suspension by means of wash centrifugation. According to the protocol of Rubinstein et al the DMSO concentration in the volume-reduced ready-to-freeze endvolume is 10% v/v. Before freezing and after thawing DMSO is potentially cytotoxic. The cytotoxicity is dependent on the DMSO concentration, the time of exposure and the temperature of the stem cell suspension during the time of exposure to DMSO (9, 10, 11, 12, 13, 14, 15, 16).

Therefore the stem cell suspension has to be kept cool at 2 °C whilst CryoSure-DEX40 is added to and removed from the stem cell suspension before freezing and after thawing, respectively. The addition of DMSO and the washing process of the stem cell suspension have to be started immediately before freezing and after thawing, respectively. In case of adequate cooling and immediate handling of the stem cell suspension during DMSO exposure in the unfrozen state, no relevant adverse effects on the cells are observed at endvolume concentrations of DMSO between 5 and 10% (15, 16, 17, 18).

Since DMSO is a strong aprotic solvent, special care has to be taken to only use DMSO-compatible materials for withdrawal of the DMSO from the vial and during transition of the DMSO to the target suspension and to minimize the contact time of DMSO with such materials. This applies even more the higher the concentration of DMSO. All processes related to the application and elimination of CryoSure-DEX40 have to be validated by the user.

The cell counting procedures carried out before freezing respectively after thawing are to be conducted by the user applying validated methods.

Storage:

CryoSure-DEX40 has a shelf life of 24 months after production if stored at 2 °C to 8 °C.

Transport:

CryoSure-DEX40 is shipped in Styrofoam boxes with ice packs to maintain a low shipping temperature. Based on stability data at elevated temperatures a transient exposure of the product to elevated temperatures during transportation should not exceed 30°C for 2 weeks. Consequently, the transient violation of the nominal storage temperature for a limited period of time, such as

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EN during transportation, does not have a negative effect on the product quality of CryoSure-DEX40. Cold chain transportation as well as temperature monitoring during transportation are therefore not necessary.

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
I Addition of CryoSure-DEX40 to the stem cell suspension

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EN I Addition of CryoSure-DEX40 to the stem cell suspension

CryoSure-DEX40 is added to the volume-reduced cord blood suspension instantaneously as the last step before initiating the freezing process.

1. Calculation of the composition of the cryoprotective solution

The amount of CryoSure-DEX40 to be added to the stem cell suspension has to be chosen in a way so that the envisaged endvolume concentration of DMSO is met. 8 ml of CryoSure-DEX40 contain 4 ml of DMSO (± 4.4 g DMSO). The density of DMSO is 1.1 g/cm^3 under normal conditions. In accordance with the protocol of Rubinstein et al 5 ml of CryoSure-DEX40 are to be added to 20 ml of volume-reduced cord blood. Thus the added 2.5 ml of DMSO result in an endvolume concentration of DMSO of 10% within the ready-to-freeze suspension (1).

2. Withdrawal of CryoSure-DEX40 from the vial and preparation of the cryoprotective solution

Before addition to the stem cell suspension CryoSure-DEX40 is to be cooled to 2°C . The amount necessary to reach the envisaged concentration of CryoSure-DEX40 in the endvolume is to be taken volumetrically from the vial.

3. Addition of the cryoprotective solution to the suspension of hematopoietic stem cells

Before addition of CryoSure-DEX40 the stem cell suspension is placed on an ice bed and cooled to 2°C . Thereafter CryoSure-DEX40, which has also been cooled to 2°C , is added volumetrically to the stem cell suspension at a constant velocity within a period of 15 minutes until the designated endvolume has been reached. Preferably a calibrated syringe pump is to be used for the addition of the cryoprotective solution (1). The consistently decelerated addition of CryoSure-DEX40 containing DMSO provides for the osmotic tolerance of the hyperosmolar DMSO and the cells in the target suspension. During the addition process, the target suspension is continuously and consistently mixed in order to assure a consistent dispensation of the DMSO conveyed within the target suspension (1).

4. Beginning of freezing process

Immediately after having added the complete amount of CryoSure-DEX40 to the target suspension, the freezing process has to be started. Until the beginning of freezing, the temperature of the stem cell suspension ready-to-freeze has to be constantly kept at 2°C . Various standard freezing procedures can be applied as specified in literature. A freezing rate of 1°C/minute until the final freeze-store temperature has been reached is considered applicable for successfully freezing hematopoietic stem cells using DMSO as a cryoprotectant (16).

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
EN II Withdrawal of cryoprotectant from the stem cell suspension after thawing

Immediately after the completed thawing process the cryoprotective solution must be washed out of the stem cell suspension. The washing process is executed in several washing steps consisting of centrifugation, withdrawal of liquid supernatant and resuspension of the cells using an appropriate wash solution. The stem cell suspension has to be kept cool at 2 °C during the washing process until the quantitative elimination of the cryoprotectant from the suspension has been concluded. Consequently the washing process has to be performed by employing a refrigerated centrifuge. The user is responsible for validating the washing process. Several methods are available for performing the washing process within a closed system (19-25).

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Sources/Quellen/Bibliographie

- 1) Rubinstein, P.; Dobrila, L.; Rosenfield, R.E.; Adamson, J.W.; Migliaccio, G.; Migliaccio, A.R.; Taylor, P.E.; Stevens, C.E. : "Processing and cryopreservation of placental / umbilical cord blood for unrelated bone marrow reconstitution" . Proc. Natl. Acad. Sci. USA 92 , 10119-10122 (1995)
- 2) Farrant 1969, Nature, 11755
- 3) Mazur 1970, Science 168, 939-949
- 4) Leibo 1977, Ciba Foundation Symposium No.52, Amsterdam S. 69-96
- 5) Merryman et al 1977 Cryobiology 14, 287-302).
- 6) Whittingham 1981, Gustav-Fischer Verlag, Stuttgart, S. 21-22
- 7) Lovelock 1953, Biochim. Biophys. Acta 10, 414-426
- 8) Merryman 1974, An. Rev. Of Biophys. And Bioengineering, Band 3 Paolo Alto, S.341-363
- 9) Goris A: '[Test of the toxicity of dimethyl sulfoxide (D.M.S.O.) on carrot tissue cultured in vitro]', Ann Pharm Fr. 1966 Dec;24(12):781-4
- 10) Basch, H., and Gadebusch, H.H. In vitro antimicrobial activity of dimethyl sulfoxide. Appl. Microbiol. 16: 1953-1954 (1968).
- 11) Chang CY, Simon E: 'The effect of dimethyl sulfoxide (DMSO) on cellular systems,' Proc Soc Exp Biol Med. 1968 May;128(1):60-6
- 12) Da Violante G, Zerrouk N, Richard I, Provot G, Chaumeil JC, Arnaud P: 'Evaluation of the cytotoxicity effect of dimethyl sulfoxide (DMSO) on Caco2/TC7 colon tumor cell cultures,' Bilo. Pharm. Bull. 25(12) 1600-1603 (2002)
- 13) Fahy GM: 'The relevance of cryoprotectant "toxicity" to cryobiology,' Cryobiology 1986 Feb;23(1):1-13
- 14) Gurtuvenko AA, Anwar J: 'Modulating the structure and properties of cell membranes: the molecular mechanism of action of dimethyl sulfoxide,' J Phys Chem B. 2007 Sep 6;111(35):10453-60
- 15) Yang H, Zhao H, Acker JP, Liu JZ, Akabutu J, McGann LE: 'Effect of dimethyl sulfoxide on post-thaw viability assessment of CD45+ and CD34+ cells of umbilical cord blood and mobilized peripheral blood,' Cryobiology 2005 Oct;51(2):165-75
- 16) Hunt CJ, Armitage SE, Pegg DE: 'Cryopreservation of umbilical cord blood: 2. Tolerance of CD34(+) cells to multimolar dimethyl sulphoxide and the effect of cooling rate on recovery after freezing and thawing,' Cryobiology. 2003 Feb;46(1):76-87
- 17) Rowley S.D., Anderson G.L.: 'Effect of DMSO exposure without cryopreservation on hematopoietic progenitor cells' Bone Marrow Transplant.11, 389-393 (1993)
- 18) Branch, D.R.; Calderwood, S.; Cecutti, M.A.; Herst, R.; Solh, H. : "Hematopoietic progenitor cells are resistant to dimethyl sulfoxide toxicity" . Transfusion 34, Nr.10 , 887-890 (1994)
- 19) Calmels B, Houzé P, Hengesse JC, Ducrot T, Malenfant C, Chabannon C.: , Preclinical evaluation of an automated closed fluid management device: Cytomate, for washing out DMSO from hematopoietic stem cell grafts after thawing,' Bone Marrow Transplant. 2003 May;31(9):823-8.
- 20) Rodríguez L, Velasco B, García J, Martín-Henao GA.: 'Evaluation of an automated cell processing device to reduce the dimethyl sulfoxide from hematopoietic grafts after thawing,' Transfusion. 2005 Aug;45(8):1391-7.
- 21) Foïs E, Desmartin M, Benhamida S, Xavier F, Vanneaux V, Rea D, Femand JP, Arnulf B, Mounier N, Ertault M, Lotz JP, Galicier L, Raffoux E, Benbunan M, Marolleau JP, Larghero J.: 'Recovery, viability and clinical toxicity of thawed and washed haematopoietic progenitor cells: analysis of 952 autologous peripheral blood stem cell transplantations,' Bone Marrow Transplant. 2007 Nov;40(9):831-5. Epub 2007 Aug 27.
- 22) Laroche V, McKenna DH, Moroff G, Schierman T, Kadidlo D, McCullough J.: , Cell loss and recovery in umbilical cord blood processing: a comparison of postthaw and postwash samples,' Transfusion. 2005 Dec;45(12):1909-16.
- 23) Lemarie C, Calmels B, Malenfant C, Arneodo V, Blaise D, Viret F, Bouabdallah R, Ladaïque P, Viens P, Chabannon C.: 'Clinical experience with the delivery of thawed and washed autologous blood cells, with an automated closed fluid management device: CytoMate,' Transfusion. 2005 May;45(5):737-42.
- 24) Nagamura-Inoue T, Shioya M, Sugo M, Cui Y, Takahashi A, Tomita S, Zheng Y, Takada K, Kodo H, Asano S, Takahashi TA.: 'Wash-out of DMSO does not improve the speed of engraftment of cord blood transplantation: follow-up of 46 adult patients with units shipped from a single cord blood bank,' Transfusion. 2003 Sep;43(9):1285-95.
- 25) Perotti CG, Del Fante C, Viarengo G, Papa P, Rocchi L, Bergamaschi P, Bellotti L, Marchesi A, Salvaneschi L.: 'A new automated cell washer device for thawed cord blood units,' Transfusion. 2004 Jun;44(6):900-6.

L-09-041-23-2-07-email (last revision: 05/2020)