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Instructions for use

(EN) CryoSure-DMSO is employed in the course of longterm storage and viability determination of cell and tissue grafts.

DMSO (Dimethyl Sulfoxide) is an intracellular cryoprotective agent which penetrates the cell membrane.

It reduces the osmotic charge of the cells during freezing and thawing and antagonizes the osmotic shock. Furthermore it protects the cells from dehydration and shrinking during freezing and prevents the formation of ice crystals.

CryoSure-DMSO is to be applied as a cryoprotectant for all cells and tissues for which the use of DMSO as a freezing agent is established.

Before transplantation CryoSure-DMSO must be eliminated from the transplant.

The procedures for viability assessment of the cell or tissue graft conducted before freezing respectively after thawing are to be validated by the user.

Each type of cell and tissue constitutes a singular system with specific characteristics regarding the behaviour during freezing and thawing. The freezing process as well as the addition and the removal of DMSO have to be validated by the user for each type of cell and each type of tissue. Amongst other factors for each cell type and each type of tissue the following parameters have to be defined:

- Mode of DMSO-addition to the transplant
- Adequate concentration of DMSO in the ready-to-freeze transplant
- Time of exposure of DMSO to the transplant and the conditions of exposure before freezing
- Freezing rate
- Storage conditions of the transplant in the frozen state
- Thawing mode and the thawing rate
- Time and conditions of exposure to DMSO of the thawed transplant before removal of the DMSO
- Mode of elimination of DMSO from the transplant before transplantation

CryoSure-DMSO may only be applied by qualified personnel who are accustomed to the use of DMSO as a cryoprotectant. Principally the application of CryoSure-DMSO consists of the following steps:

1. Withdrawal of DMSO from the sterile container
2. Addition of DMSO to the cell suspension respectively the tissue
3. Freezing of the transplant
4. Thawing of the transplant
5. Elimination of DMSO from the thawed transplant before transplantation

1. Withdrawal of DMSO from the sterile container:

The necessary amount of DMSO for freezing is taken volumetrically from the sterile container. The opening of the sterile container and the withdrawal of DMSO thereof shall be executed under conditions as to minimize the risk of contamination. Since DMSO is a strong aprotic solvent which solves hydrophobic as well as hydrophilic substances it is imperative to take care only to use

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materials for withdrawal and transition of DMSO which are DMSO-compatible and to keep the contact time with these materials as short as possible.

2. Addition of DMSO to the cell suspension respectively the tissue:

During addition of DMSO the following factors have to be considered:

Release of heat:

During addition of DMSO to water heat is released. This heat must be dissipated before contact with the cells or tissue in order to avoid negative impact of the heat released. This can be achieved by preparing a cryoprotective solution of DMSO and an appropriate vector solution. The release of heat can be dissipated via an ice bed. Such vector solutions can be for example autologous plasma, plasma-like solutions, cell culture media or such vector solutions which contain other cryoprotectants. The release of heat during DMSO-addition may lead to the reduction of the volume of the aqueous target solution.

Osmotic charge of cells:

DMSO has a high osmolality. Therefore the addition of DMSO to cells or tissue leads to outstream of water from the cell and therefore to the shrinking of the cells until the intracellular and extracellular concentrations of DMSO have equilibrated. The optimal mode of addition is specific for each type of cell and each type of tissue and has to be validated by the user. It has to be taken care that the cells do not shrink over a critical point during the DMSO-addition. Preferably a calibrated syringe pump is to be used for the DMSO-addition.

Cytotoxicity of DMSO:

Prior to freezing respectively after thawing DMSO is potentially cytotoxic. The cytotoxicity of DMSO is dependent on the DMSO-concentration, the time of exposure and the temperature of the cell suspension or the tissue during the exposure. The conditions of DMSO-addition and the conditions of DMSO-exposure to the transplant (e.g. time and temperature) are to be validated and defined by the user for each cell type and type of tissue. The concentration of DMSO in a ready-to-freeze cell suspension is typically between 5 and 10 % v/v. The DMSO-concentration in the solution in which the tissue transplant is placed is typically higher than in the cell suspension, e.g. up to 20 % v/v. The temperature of the ready-to-freeze cell suspension respectively tissue during DMSO-addition should be kept at 2 - 4 °C. The time of DMSO-exposure to the cells respectively the tissue in the unfrozen state should be kept as short as possible. In mononuclear blood cells for example the exposure time should not exceed 15 minutes. In case of tissues a longer exposure is eventually necessary in order to achieve the necessary DMSO concentration throughout the tissue. In such cases the exposure time should not exceed one hour. In case the transplant is processed at higher temperatures than indicated above it has to be taken into

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consideration that the cytotoxicity of DMSO increases with the temperature.

3. Freezing of the transplant:

The optimum freezing rate is dependent on the cell type and type of tissue and has to be validated and defined by the user. The freezing is to be executed according to a standardized procedure, preferably with a controlled rate freezer.

4. Thawing of the transplant:

The thawing mode and the thawing rate are specific for each cell type and type of tissue and have to be validated and defined by the user.

5. Elimination of DMSO from the thawed transplant before transplantation:

After having finished the thawing process DMSO must be eliminated from the transplant. The elimination process is executed by means of washing and dilution steps employing adequate wash solutions and equipment. The elimination mode is specific to each cell type and type of tissue and has to be validated and defined by the user. Generally the following factors have to be considered:

Osmotic charge of cells:

Contrary to the addition process each washing step initially results in an elevated intracellular DMSO-concentration compared to the extracellular concentration. This leads to an inflow of water into the cell and consequently the cell volume increases until the intracellular and extracellular DMSO-concentrations have equilibrated. The mode of DMSO-elimination is to be chosen in a way so that the cell volume does not increase over a critical point. The optimal elimination mode is specific for each cell type and type of tissue and is to be validated and defined by the user.

Cytotoxicity of DMSO:

Before freezing and after thawing DMSO is potentially cytotoxic. The cytotoxicity is dependent on the DMSO-concentration, the time of exposure and the temperature of the cell suspension or tissue during the DMSO-exposure. Analogously to the conditions before freezing the time of exposure of the transplant to DMSO in the non-frozen state should be kept as short as possible, and the temperature during the elimination of DMSO should be kept between 2 °C and 4 °C. The conditions of DMSO-elimination and the conditions of DMSO-exposure to the thawed transplant (e.g. time and temperature) are to be validated and defined by the user.

6. General directions to the use of CryoSure-DMSO:

• Penetration of DMSO via the skin:

DMSO has the characteristic to penetrate via the human skin and then to act systemically. Furthermore certain alien substances can be carried along across the skin barrier and eventually also act systemically. Therefore it is necessary to wear protective gloves, which are non-permeable to DMSO, during handling of CryoSure-DMSO and to take care that

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the DMSO does not get in touch with substances that could be carried via the skin barrier by DMSO and that could consequently have negative effects.

- **Utilization of DMSO-compatible materials:**

Since DMSO is a strong aprotic solvent, it is necessary to take care that during withdrawal from the container and during transition DMSO only comes in contact with materials which are compatible with DMSO and to keep the contact time with such materials as short as possible. The higher the DMSO-concentration materials get in contact with, the shorter the contact time with such materials should be kept.

- **Melting point of DMSO:**

Pure DMSO has its melting point at 18.3 °C. At lower temperatures DMSO solidifies. The change of the aggregate state within the closed sterile container does not affect the product quality. However CryoSure-DMSO should only be applied as a cryoprotectant in the liquid state. In order to regain the liquid state of DMSO CryoSure-DMSO should be kept at 20 - 30 °C until the liquid state is fully reconstituted prior to application.

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