

Product Description

The dNTP Mix is an equimolar mixture of high purity dATP, dTTP, dCTP and dGTP, which can be directly used for PCR, reverse transcription (RT) and one-step RT-PCR experiment.

Components

Components	P031-01	P031-02
dNTP Mix (10 mM each)	1 ml	5 × 1 ml

▲ The above components are sodium saline solution of dNTP, pH 7.0 - 9.0.

Storage

Store at -30 ~ -15°C and transport at ≤0°C.

Applications

It is applicable for PCR, reverse transcription (RT) and one-step RT-PCR experiment.

Notes

For research use only. Not for use in diagnostic procedures.

Experiment Process

This mix can be used directly without dilution.

For a 50 µl PCR system, it is recommended to add 1 µl of P031-01/02.

For a 20 µl RT system, it is recommended to add 1 µl of P031-01/02.