

Certificate of Quality

PRODUCT: 10X PBS Buffer with Tween-20

Cat. No.: BR540-1L

Lot No.:

COMPONENTS: 1.370M NaCl
0.100M Na₂HPO₄
0.027M KCl
0.018M KH₂PO₄
0.5% Tween 20
pH 7.4

DESCRIPTION: Results: Pass

STORAGE: Store at room temperature or at 4°C for longer storage.

CAUTION: Avoid contact with eyes, skin and clothing.
Wash all areas of contact with copious amounts of water.
For laboratory or further manufacturing use only.